

Installation Instructions

Transformers (Field Wired)

Part #ACA110, ACA112

Description:

Apollo Dental Products' Transformers are highly recommended in areas of inadequate electrical line voltage. Buck-Boost Transformers provide an economical and convenient means of either increasing or decreasing available line voltage for 115 or 230 volt electrical circuits. The ADP Transformer should be installed by a qualified electrician.

DIMENSIONS:		
Model	ACA110	ACA112
Width	7 7/8"	4 1/2"
Depth	5 1/2"	4 1/2"
Height	9 5/8"	7 3/4"
Weight	25 lbs.	20 lbs.
SPECIFICATIONS:		
Voltage	115/230	230
KVA	7.2	4.7
Boost Volts	96-115 / 208-230	208-230
Connection	Field Wired	Field Wired

Warranty Information: 90 Days

All ADP accessories are thoroughly inspected and tested in accordance with rigid specifications and standards. Our accessories are guaranteed against any defective material and workmanship from the date of shipment; provided that the installation, operation, and maintenance is done in accordance with ADP procedures as outlined in our Installation and Maintenance Guides. Warranty cards must be returned to ADP within ten days of installation to effect warranty. No other warranties or guarantees, expressed or implied are made.

ADP's obligation under the warranty is to provide parts for the repair or, at its option, to provide the replacement product (excluding labor). All special, incidental and/or consequential damages are excluded. We will not issue credit for accessories without first attempting to correct the problem in the field. Written notice of breach of warranty must be given to ADP within the warranty period. The warranty does not cover damage resulting from improper installation or maintenance, accident or misuse. The warranty does not cover damage resulting from the use of cleaning, disinfecting or sterilizing chemicals and processes. Failure to follow instructions provided in ADP's Installation and Maintenance Guides may void the warranty.

Symbols:

X - # = Wire Code On Transformer

← = Wire Nut Together

—●— = Continue On To Compressor Or Device And Ties To Wire On Transformer

Note: Ground Wire Should Attach To Housing Of Transformer.

